

USP28 Rabbit mAb


ABZ1897

Product Overview

Name	USP28 Rabbit mAb
Accession(Primary)	Q96RU2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:100-1:200

Target information (Q96RU2)

Gene Name	USP28
About Gene Symbol	



Western blot analysis of USP28 in A549 lysates
using USP28 antibody.