

RAP80 Rabbit mAb


ABZ1899

Product Overview


Name	RAP80 Rabbit mAb
Accession(Primary)	Q96RL1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 80 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (Q96RL1)


Gene Name	UIMC1
About Gene Symbol	



Immunocytochemistry analysis of RAP80 (green) in HeLa using RAP80 antibody, and DAPI (blue).



Western blot analysis of RAP80 in HeLa, A549, HL-60, U2OS lysates using RAP80 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using RAP80 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.