

TRPM7 Rabbit mAb


ABZ1900

Product Overview

Name	TRPM7 Rabbit mAb
Accession(Primary)	Q96QT4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 213 kDa; Observed MW: 213 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (Q96QT4)

Gene Name	TRPM7
About Gene Symbol	TRPM7 (Transient Receptor Potential Cation Channel Subfamily M Member 7) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of TRPM7 in HepG2, JAR, lysates using TRPM7 antibody.



Western blot analysis of TRPM7 in HeLa lysates using TRPM7 antibody.