

SMG1 Rabbit mAb

ABZ1902


Product Overview

Name	SMG1 Rabbit mAb
Accession(Primary)	Q96Q15
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 411 kDa; Observed MW: 411 kDa
Recommended Dilution	WB-1:1000-1:5000


Target information (Q96Q15)

Gene Name SMG1

About Gene Symbol



Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Smg1 in Saos2 lysates
using SMG1 antibody.