

QK1 Rabbit mAb


ABZ1904

Product Overview


Name	QK1 Rabbit mAb
Accession(Primary)	Q96PU8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; IP-1:10-1:100

Target information (Q96PU8)

Gene Name	QKI
About Gene Symbol	QKI (QKI, KH Domain Containing RNA Binding) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of QK1 in 3T3, rat brain lysates using QK1 antibody.



Western blot analysis of QK1 in K562 lysates using QK1 antibody.