

Importin 9 Rabbit mAb


ABZ1906

Product Overview


Name	Importin 9 Rabbit mAb
Accession(Primary)	Q96P70
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:200-1:1000; ICC/IF-1:50-1:100

Target information (Q96P70)

Gene Name	IPO9
About Gene Symbol	IPO9 (Importin 9) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of Importin 9 (green) in A549 using Importin 9 antibody, and DAPI(blue).



Western blot analysis of Importin 9/RANBP9 in Jurkat, HeLa lysates using Importin 9/RANBP9 antibody.