


NLRP3 Rabbit mAb

ABZ1907

Product Overview

Name	NLRP3 Rabbit mAb
Accession(Primary)	Q96P20
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 118 kDa; Observed MW: 118 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96P20)

Gene Name NLRP3**About Gene Symbol**

Western blot analysis of NLRP3 in THP-1 lysates
using NLRP3 antibody.

