

PHD1/prolyl hydroxylase Rabbit mAb


ABZ1908

Product Overview


Name	PHD1/prolyl hydroxylase Rabbit mAb
Accession(Primary)	Q96KS0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (Q96KS0)

Gene Name	EGLN2
About Gene Symbol	



Western blot analysis of PHD1/prolyl hydroxylase in rat heart lysates using PHD1 antibody.



Western blot analysis of PHD1/prolyl hydroxylase in HeLa, A549, HL-60, U2OS, U87-MG lysates using PHD1/prolyl hydroxylase antibody.