

PBK Rabbit mAb


ABZ1909

Product Overview


Name	PBK Rabbit mAb
Accession(Primary)	Q96KB5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 36 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IP-1:20-1:50

Target information (Q96KB5)

Gene Name	PBK
About Gene Symbol	PBK (PDZ Binding Kinase) is a serine/threonine protein kinase involved in MAPK/ERK signaling and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded mouse colon using PBK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PBK in K562, HeLa lysates using PBK antibody.