


MS4A14 Rabbit mAb

ABZ1912

Product Overview

Name	MS4A14 Rabbit mAb
Accession(Primary)	Q96JA4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q96JA4)

Gene Name MS4A14**About Gene Symbol**

Western blot analysis of MS4A14 in Jurkat, rat Brain, C6 lysates using MS4A14 antibody.

