


ITCH/AIP4 Rabbit mAb

ABZ1913

Product Overview

Name	ITCH/AIP4 Rabbit mAb
Accession(Primary)	Q96J02
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96J02)

Gene Name ITCH**About Gene Symbol**

Western blot analysis of ITCH/AIP4 in A549, HL-60 lysates using ITCH antibody.

