


HPDL Rabbit mAb

ABZ1914

Product Overview

Name	HPDL Rabbit mAb
Accession(Primary)	Q96IR7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q96IR7)

Gene Name HPDL**About Gene Symbol**

Western blot analysis of HPDL in HeLa, Caco2, A549 lysates using HPDL antibody.

