

Pellino 1 Rabbit mAb

ABZ1920


Product Overview

Name	Pellino 1 Rabbit mAb
Accession(Primary)	Q96FA3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20; IF-1:50-1:200

Target information (Q96FA3)

Gene Name PELI1

About Gene Symbol



Western blot analysis of Pellino 1 in rat Brain, Jurkat lysates using Pellino 1 antibody.

