

Nogo B Receptor Rabbit mAb

ABZ1921


Product Overview

Name	Nogo B Receptor Rabbit mAb
Accession(Primary)	Q96E22
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96E22)

Gene Name NUS1

About Gene Symbol



Western blot analysis of Nogo B receptor in K562, 3T3, HeLa lysates using Nogo B Receptor antibody.

