

SLAMF6 Rabbit mAb


ABZ1922

Product Overview

Name	SLAMF6 Rabbit mAb
Accession(Primary)	Q96DU3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 37 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (Q96DU3)

Gene Name	SLAMF6
About Gene Symbol	



Western blot analysis of SLAMF6 in Ramos lysates using SLAMF6 antibody.