

Optineurin Rabbit mAb


ABZ1924

Product Overview


Name	Optineurin Rabbit mAb
Accession(Primary)	Q96CV9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 66 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:500-1:1000; FC-1:20-1:50; IP-1:20-1:100

Target information (Q96CV9)

Gene Name	OPTN
About Gene Symbol	OPTN (Optineurin) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Optineurin in HeLa, A549, HL-60, U2OS, U87-MG lysates using Optineurin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Optineurin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.