

Rab7B Rabbit mAb


ABZ1926

Product Overview

Name	Rab7B Rabbit mAb
Accession(Primary)	Q96AH8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; FC-1:100-1:200

Target information (Q96AH8)

Gene Name	RAB7B
About Gene Symbol	



Western blot analysis of RAB7B in HeLa lysates
using Rab7B antibody.