

## FUBP1/FBP Rabbit mAb


ABZ1927

### Product Overview


|                             |                                                                                           |
|-----------------------------|-------------------------------------------------------------------------------------------|
| <b>Name</b>                 | FUBP1/FBP Rabbit mAb                                                                      |
| <b>Accession(Primary)</b>   | Q96AE4                                                                                    |
| <b>Host Species</b>         | Rabbit                                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                                |
| <b>Isotype</b>              | IgG                                                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                         |
| <b>Calculated MW</b>        | Calculated MW: 68 kDa; Observed MW: 74 kDa                                                |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:200-1:500; FC-1:200-1:500; IP-1:20-1:50 |

### Target information (Q96AE4)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | FUBP1 |
| <b>About Gene Symbol</b> |       |




Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).



Western blot analysis of FUBP1/GBP in mouse spleen lysates using FUBP1/GBP antibody.



Western blot analysis of FUBP1/GBP in K562, rat Brain, HeLa lysates using FUBP1/GBP antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.