

FUBP1/FBP Rabbit mAb


ABZ1928

Product Overview


Name	FUBP1/FBP Rabbit mAb
Accession(Primary)	Q96AE4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (Q96AE4)


Gene Name	FUBP1
About Gene Symbol	




Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).



Western blot analysis of FUBP1/FPB in mouse brain, mouse spleen lysates using FUBP1/FPB antibody.



Western blot analysis of FUBP1/GBP in Jurkat, C6 lysates using FUBP1/GBP antibody



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.