

Adipose Triglyceride Lipase Rabbit mAb


ABZ1929

Product Overview


Name	Adipose Triglyceride Lipase Rabbit mAb
Accession(Primary)	Q96AD5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:200-1:500; IP-1:20-1:50

Target information (Q96AD5)

Gene Name	PNPLA2
About Gene Symbol	



Western blot analysis of Adipose Triglyceride Lipase in K562, C6, 3T3, HeLa lysates using Adipose Triglyceride Lipase antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Adipose Triglyceride Lipase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.