

Adipose Triglyceride Lipase Rabbit mAb

ABZ1930


Product Overview


Name	Adipose Triglyceride Lipase Rabbit mAb
Accession(Primary)	Q96AD5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q96AD5)

Gene Name PNPLA2

About Gene Symbol





Western blot analysis of Adipose Triglyceride Lipase in K562 lysates using Adipose Triglyceride Lipase antibody.