

Wnt2b Rabbit mAb


ABZ1933

Product Overview


Name	Wnt2b Rabbit mAb
Accession(Primary)	Q93097
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:100

Target information (Q93097)

Gene Name	WNT2B
About Gene Symbol	



Immunocytochemistry analysis of Wnt2b (green) in 293t using Wnt2b antibody, and DAPI(blue).



Western blot analysis of Wnt2b in K562, C6, 3T3, HeLa lysates using Wnt2b antibody.