

Cullin 5/CUL-5 Rabbit mAb


ABZ1934

Product Overview

Name	Cullin 5/CUL-5 Rabbit mAb
Accession(Primary)	Q93034
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q93034)

Gene Name	CUL5
About Gene Symbol	CUL5 (Cullin 5) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cullin 5/CUL5 in C6, 3T3, HeLa lysates using Cullin 5 antibody.