

## USP7 Rabbit mAb


ABZ1935

### Product Overview


|                             |                                                                       |
|-----------------------------|-----------------------------------------------------------------------|
| <b>Name</b>                 | USP7 Rabbit mAb                                                       |
| <b>Accession(Primary)</b>   | Q93009                                                                |
| <b>Host Species</b>         | Rabbit                                                                |
| <b>Clonality</b>            | Monoclonal                                                            |
| <b>Isotype</b>              | IgG                                                                   |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                     |
| <b>Calculated MW</b>        | Calculated MW: 128 kDa; Observed MW: 128 kDa                          |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC:1:50-1:100;<br>ICC/IF-1:100-1:200; FC-1:20-1:50 |

### Target information (Q93009)

|                          |                                                                                                                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | USP7                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | USP7 (Ubiquitin Specific Peptidase 7) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of USP7 in C6 lysates using USP7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HAUSP/USP7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.