

USP13 Rabbit mAb


ABZ1936

Product Overview


Name	USP13 Rabbit mAb
Accession(Primary)	Q92995
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:1000-1:5000

Target information (Q92995)

Gene Name	USP13
About Gene Symbol	USP13 (Ubiquitin Specific Peptidase 13) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of USP13 (green) in Hela using USP13 antibody, and DAPI (blue).



Western blot analysis of USP13 in U2OS lysates using USP13 antibody.