

DLX4 Rabbit mAb


ABZ1937

Product Overview



Name	DLX4 Rabbit mAb
Accession(Primary)	Q92988
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; IHC-1:50-1:100

Target information (Q92988)

Gene Name	DLX4
About Gene Symbol	



Immunocytochemistry analysis of DLX4 (green) in K562 using DLX4 antibody, and DAPI (blue).



Western blot analysis of Dlx4 in K562, rat Brain, Hela lysates using Dlx4 antibody.