

IRF7 Rabbit mAb


ABZ1938

Product Overview


Name	IRF7 Rabbit mAb
Accession(Primary)	Q92985
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; FC-1:50-1:100; IP-1:10-1:100

Target information (Q92985)


Gene Name	IRF7
About Gene Symbol	



Immunocytochemistry analysis of IRF7 (green) in 293 using IRF7 antibody, and DAPI (blue).



Western blot analysis of IRF7 in rat Brain lysates using IRF7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using IRF7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.