

Bad Rabbit mAb


ABZ1939

Product Overview


Name	Bad Rabbit mAb
Accession(Primary)	Q92934
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q92934)


Gene Name	BAD
About Gene Symbol	



Western blot analysis of Bad in rat Brain, C6, 3T3, Hela lysates using Bad antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Bad antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Bad (green) in Hela using Bad antibody, and DAPI (blue)

