

UFD1L Rabbit mAb


ABZ1942

Product Overview

Name	UFD1L Rabbit mAb
Accession(Primary)	Q92890
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (Q92890)

Gene Name	UFD1
About Gene Symbol	



Western blot analysis of UFD1L in HeLa lysates
using UFD1 antibody.