


Protein Jumonji Rabbit mAb

ABZ1948

Product Overview

Name	Protein Jumonji Rabbit mAb
Accession(Primary)	Q92833
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 139 kDa; Observed MW: 139 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100

Target information (Q92833)

Gene Name JARID2**About Gene Symbol**

Western blot analysis of Jarid2 in K562, C6 lysates using Protein Jumonji antibody.

