

Golgin 97 Rabbit mAb

ABZ1950


Product Overview

Name	Golgin 97 Rabbit mAb
Accession(Primary)	Q92805
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q92805)

Gene Name GOLGA1

About Gene Symbol



Western blot analysis of Golgin 97 in HeLa lysates using GOLGA1 antibody.

