

PROX1 Rabbit mAb

ABZ1951


Product Overview

Name	PROX1 Rabbit mAb
Accession(Primary)	Q92786
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q92786)

Gene Name PROX1

About Gene Symbol



Western blot analysis of PROX1 in C6 lysates using PROX1 antibody.

