

HDAC2 Rabbit mAb


ABZ1953

Product Overview


Name	HDAC2 Rabbit mAb
Accession(Primary)	Q92769
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:500-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (Q92769)


Gene Name	HDAC2
About Gene Symbol	HDAC2 (Histone Deacetylase 2) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC2 in K562, C6, 3T3 lysates using HDAC2 antibody.



Immunocytochemistry analysis of HDAC2(green) in HeLa using HDAC2 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using HDAC2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.