

AP2 gamma/TFAP2C Rabbit mAb


ABZ1954

Product Overview


Name	AP2 gamma/TFAP2C Rabbit mAb
Accession(Primary)	Q92754
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:500-1:1000; HC-P:1:200-1:2000; ICC/IF-1:200-1:500; FC-1:20-1:50; IP-1:20-1:50; IHC-1:50-1:100

Target information (Q92754)

Gene Name	TFAP2C
About Gene Symbol	TFAP2C (Transcription Factor AP-2 Gamma) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TFAP2C/AP2 gamma in Jurkat, HeLa lysates using Transcription Factor AP 2 gamma antibody.



Immunocytochemistry analysis of TFAP2C/AP2 gamma (green) in MCF-7 using TFAP2C/AP2 gamma antibody, and DAPI (blue)