

## AP2 gamma/TFAP2C Rabbit mAb


ABZ1955

### Product Overview


|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | AP2 gamma/TFAP2C Rabbit mAb                                               |
| <b>Accession(Primary)</b>   | Q92754                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 49 kDa; Observed MW: 49 kDa                                |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:200-1:500; FC-1:50-1:100 |

### Target information (Q92754)


|                          |                                                                                                                                                                                                                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TFAP2C                                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | TFAP2C (Transcription Factor AP-2 Gamma) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Immunocytochemistry analysis of Transcription Factor AP2 gamma (green) in HeLa using Transcription Factor AP2 gamma antibody, and DAPI (blue).



Western blot analysis of AP2 gamma/TFAP2C in HeLa, A549 lysates using AP2 gamma/TFAP2C antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using AP2 gamma/TFAP2C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.