

TFG Rabbit mAb


ABZ1956

Product Overview


Name	TFG Rabbit mAb
Accession(Primary)	Q92734
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500

Target information (Q92734)


Gene Name	TFG
About Gene Symbol	TFG (Trafficking From ER To Golgi Regulator) is a protein kinase and signaling regulator involved in vesicular trafficking. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TRK fused gene in C6, 3T3 lysates using TFG antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TRK fused gene antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TRK fused gene (green) in Hela using TRK fused gene antibody, and DAPI(blue)