

Nectin 2 Rabbit mAb


ABZ1957

Product Overview


Name	Nectin 2 Rabbit mAb
Accession(Primary)	Q92692
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 70-80 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:20; IP-1:20

Target information (Q92692)

Gene Name	NECTIN2
About Gene Symbol	



Western blot analysis of Nectin 2 in HeLa, A549, HL-60, U251, U87-MG lysates using Nectin 2 antibody.



Immunocytochemistry analysis of Nectin 2 (green) in HeLa using Nectin 2 antibody, and DAPI (blue)