

PHAPI2 Rabbit mAb


ABZ1958

Product Overview



Name	PHAPI2 Rabbit mAb
Accession(Primary)	Q92688
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:20

Target information (Q92688)


Gene Name	ANP32B
About Gene Symbol	ANP32B (Acidic Nuclear Phosphoprotein 32 Family Member B) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of PHAPI2 (green) in K562 using PHAPI2 antibody, and DAPI (blue).



Western blot analysis of PHAPI2 in K562, C6 lysates using PHAPI2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer tissue using PHAPI2/APRIL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.