

Neurogranin Rabbit mAb


ABZ1959

Product Overview

Name	Neurogranin Rabbit mAb
Accession(Primary)	Q92686
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 8 kDa; Observed MW: 14-16 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; ICC/IF-1:500-1:2000; IP-1:20-1:50

Target information (Q92686)

Gene Name	NRGN
About Gene Symbol	



Western blot analysis of Neurogranin in U87-MG, PC-3 lysates using Neurogranin antibody.