


MRPS31 Rabbit mAb

ABZ1960

Product Overview

Name	MRPS31 Rabbit mAb
Accession(Primary)	Q92665
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q92665)

Gene Name MRPS31**About Gene Symbol**

Western blot analysis of MRPS31 in rat Brain lysates using MRPS31 antibody.

