

PI3 Kinase p55 gamma Rabbit mAb


ABZ1962

Product Overview


Name	PI3 Kinase p55 gamma Rabbit mAb
Accession(Primary)	Q92569
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50; IHC-1:50-1:100

Target information (Q92569)


Gene Name	PIK3R3
About Gene Symbol	




Western blot analysis of PI3 Kinase p55 gamma in HEK-293T, 3T3, C6 lysates using PI3 Kinase p55 gamma antibody.



Western blot analysis of PI 3 Kinase p55 gamma in SH-SY5Y lysates using PI3 Kinase p55 gamma antibody.



Immunohistochemistry analysis of paraffin-embedded Human astrocytoma using PI 3 Kinase p55 gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of PI3 Kinase p55 gamma in C6 using PI 3 Kinase p55 gamma antibody.