


FIG4 Rabbit mAb

ABZ1963

Product Overview

Name	FIG4 Rabbit mAb
Accession(Primary)	Q92562
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 104 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q92562)

Gene Name FIG4**About Gene Symbol**

Western blot analysis of FIG4 in HL-60, U2OS, U87-MG lysates using FIG4 antibody.

