

BAP1 Rabbit mAb


ABZ1964

Product Overview

Name	BAP1 Rabbit mAb
Accession(Primary)	Q92560
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 80 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:100; FC-1:50-1:100

Target information (Q92560)

Gene Name	BAP1
About Gene Symbol	



Western blot analysis of BAP1 in HeLa, A549, HL-60, U2OS lysates using BAP1 antibody.