


WAVE1 Rabbit mAb

ABZ1965

Product Overview

Name	WAVE1 Rabbit mAb
Accession(Primary)	Q92558
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (Q92558)

Gene Name WASF1**About Gene Symbol**

Western blot analysis of WAVE 1 in rat Brain, C6, 3T3, HeLa lysates using WAVE 1 antibody.

