

ELMO1 Rabbit mAb


ABZ1966

Product Overview


Name	ELMO1 Rabbit mAb
Accession(Primary)	Q92556
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:10-1:100

Target information (Q92556)


Gene Name	ELMO1
About Gene Symbol	



Immunocytochemistry analysis of ELMO1 (green) in Jurkat using ELMO1 antibody, and DAPI (blue).



Western blot analysis of ELMO1 in Jurkat lysates using ELMO1 antibody.



Western blot analysis of ELMO1 in rat Brain, C6,
3T3 lysates using ELMO1 antibody