

SIRT1 Rabbit mAb


ABZ1967

Product Overview

Name	SIRT1 Rabbit mAb
Accession(Primary)	Q923E4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 80 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-P-1:200-1:500; IF: 1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (Q923E4)

Gene Name	Sirt1
About Gene Symbol	



Western blot analysis of SIRT1 in K562, C6, 3T3, HeLa lysates using SIRT1 antibody.