

DGCR8 Rabbit mAb


ABZ1968

Product Overview



Name	DGCR8 Rabbit mAb
Accession(Primary)	Q8WYQ5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; FC-1:200-1:1000; IP-1:20-1:50

Target information (Q8WYQ5)

Gene Name	DGCR8
About Gene Symbol	DGCR8 (DGCR8 Microprocessor Complex Subunit) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of DGCR8 (green) in K562 using DGCR8 antibody, and DAPI (blue).



Western blot analysis of DGCR8 in K562, C6, 3T3 lysates using DGCR8 antibody.