


ATG4A Rabbit mAb

ABZ1969

Product Overview

Name	ATG4A Rabbit mAb
Accession(Primary)	Q8WYN0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q8WYN0)

Gene Name ATG4A**About Gene Symbol**

Western blot analysis of ATG4A in Jurkat lysates using ATG4A antibody.

