

BLNK Rabbit mAb

ABZ1971


Product Overview

Name	BLNK Rabbit mAb
Accession(Primary)	Q8WV28
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q8WV28)

Gene Name BLNK

About Gene Symbol



Western blot analysis of BLNK in K562, C6, 3T3, HeLa lysates using BLNK antibody.

