


PHLDA1 Rabbit mAb

ABZ1972

Product Overview

Name	PHLDA1 Rabbit mAb
Accession(Primary)	Q8WV24
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q8WV24)

Gene Name PHLDA1**About Gene Symbol**

Western blot analysis of PHLDA1 in A375 lysates using PHLDA1 antibody.

